




CITY 2 PLUS by Travel Buggy

A HEAVY DUTY AND PREMIUM FOLDING POWER CHAIR

NOW EASIER TO LIFT

3-PIECE ASSEMBLY



HILLS DONT STAND A CHANCE

The CITY 2 PLUS is our strongest and most powerful chair in our line-up. It features dual batteries for an extra long driving range, a reinforced aluminum frame capable of holding up to 395 lbs, and powerful motors that can tackle any hill.



DRIVING RANGE
15.5 MI / 25 KM



WEIGHT CAPACITY
395"



INCLINE ABILITY
12.0° / 21.3%



ADJUSTABLE SEATING FOR ULTIMATE COMFORT

Get comfortable in your CITY 2 PLUS with a large economically designed seat cushion, adjustable seating width & soft contoured arm rests. Utilize the reclining seat back to reduce back pain and extend the footrest up to 5" to help with sore knees.



RECLINING BACK
115°




WIDENING SEAT
19.5" - 23"



FOOT PLATE
ADJUSTABLE 5"

Visit your Authorised Travel Buggy Dealer:



Intelligently Designed

KEY SPECIFICATIONS



GROSS/NET WEIGHT
65 LBS / 58 LBS



MOTOR POWER
2 X 250 W BRUSHLESS



BATTERY TYPE
AIRLINE APPROVED



MAX SPEED
4.5 MPH / 7 KPH



FOLDED DIMENSIONS
33" L X 23" W X 13" H



UNFOLDED DIMENSIONS
44" L X 23" W X 40" H



PERFECT FOR
INDOOR & OUTDOOR USE




REAR WHEELS
LARGE 12"



When you choose Travel Buggy, you're not just getting a mobility solution; you're becoming part of a vibrant community of adventurers, explorers, and advocates for accessible travel. It is time to **SEE YOUR WORLD AGAIN** with Travel Buggy.



TRAVEL BUGGY WAS FOUNDED TO PROVIDE
TRAVEL-FRIENDLY POWER MOBILITY EQUIPMENT



We offer more accessories than anyone else. From essentials like cup holders & extra storage to rehab accessories like specialized seating, head rests & leg extensions. We can tailor your AEROLUX to your desire!

