




Standard Colors:



Crimson Red



Stone Grey



Desert Tan

Eagle 4-Wheel GR595



If you're looking for the ultimate outdoor driving experience, look no further than the new Golden Eagle™. In a class all its own, the Eagle's stylish yet rugged, off-road design is built for high performance and maneuverability with front and rear suspension for smoothly traversing outdoor terrain and a powerful 140 Amp controller for outstanding climbing ability. It offers a full LED lighting package, an adjustable tiller and a luxury high-back Captain's seat for maximum comfort when on the go, plus durable vibrant red color shrouds for a sporty look. With a 500 lb. weight capacity, 13" pneumatic tires, up to 32-mile operating range on optional 75 AH batteries, plus a water-resistant storage compartment that comes standard, the Golden Eagle offers everything you need to enjoy the great outdoors!

Standard Features:

- 500 lb. Weight Capacity
- Integrated LCD Console Display
- Full LED Lighting Package
- Built-in Speaker with Bluetooth® & USB Connectivity
- Lockable Water-resistant Storage Compartment
- Rearview Mirrors
- Cell Phone Holder & USB Charger
- Convenient Drink Holders
- Infinite Position Tiller
- High-back, Limited Recline Captain's Seat
- Large Front Basket
- 13" Pneumatic Tires
- 2-NF22 Batteries

**Upgrade with
75 AH Batteries for
Extended Range!**



| Item | GR595 |
|--|---------------------------------------|
| Medicare Code | Not Coded |
| Weight Capacity | 500 lbs |
| Drive Wheels | Rear |
| Maximum Speed | 8 mph |
| Operating Range ¹ | 23 miles / 32 miles |
| Ground Clearance | 5" |
| Turning Radius | 69" |
| Type Batteries | 2-NF22 / 2-75 AH |
| Color | Red, Grey & Tan |
| Freewheel Mode | Yes |
| Electronic Speed Control | Yes |
| Electro-Mechanical Brakes | Yes |
| Mechanical Hand Brake | No |
| Charger | Off Board |
| Controller | 140 Amp |
| Dimensions: | |
| Length | 56.5" |
| Width | 27" |
| Height (ground to top of back of seat) | 49" - 50" |
| Height (ground to top of back of headrest) | 57" - 62.5" |
| Ground to Top of Seat | 27" - 28" |
| Ground to Top of Deck | 7.25" |
| Deck to Top of Seat | 19.75" - 20.75" |
| Number of Seat Height Positions | 2 |
| Size of Increments | 1" |
| Front Axle to Rear Axle | 40" |
| Weight of Unit: (assembled) | |
| Base | 179 lbs |
| Seat w/Arms | 49.5 lbs |
| Batteries | 37 lbs ea (NF22) 52 lbs ea (75 AH) |
| Tires: | |
| Front ³ | Pneumatic 13" |
| Rear ³ | 13" |
| Standard Seat: | |
| Slide Seat | Captain's High-Back Standard |
| Back Height - no headrest | 23.75" |
| Back Height with headrest | 30.75" - 36.25" |
| Width x Depth | 20" x 18" |
| Color | Black |
| Warranty: | |
| Frame | 5 Years |
| Drive Train | 2 Year |
| Electronics ² | 2 Year |
| Batteries | 1 Year |

¹ Battery range at 200 lbs., and will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.

³ Inflate tires to 35 psi.



Rear view mirrors, cell phone holder and built-in speaker system with Bluetooth® & USB connectivity.



Integrated LCD console display monitors speed, odometer, fault codes & battery life.



Lockable water-resistant storage compartment comes standard and attaches to the accessory bracket.



Ultrabright LED rear lights.



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com

Find us on Facebook

www.facebook.com/goldentechnologies