

AL600 CLASSIC

THE ULTIMATE INSIDE PLATFORM LIFT FOR
SCOOTERS AND POWERCHAIRS



KEEPING IT INSIDE

The AL600 Classic is an ideal solution to keep your scooter or powerchair protected inside your vehicle during transit. Guide your mobility device on and off the platform for easy loading/unloading. At night, enjoy the new, backlit remote control with a built-in flashlight! Say "goodbye" to struggling with straps in the dark. Don't be fooled by its compact dimensions - the lift boasts a 35.5" vertical travel from the cargo area to the ground with the press of a button.

Opt for the best of both worlds of the hoist and hitch mount lifts: Discreet mobility device transportation, continued use of your vehicle's backup camera, and the ease of a platform lift all in one! Designed to do the heavy lifting for you, the AL600 Classic is the optimal choice for transporting your mobility device inside vans and SUVs.

EASY ON, EASY OFF ACCESS

AL600 CLASSIC



The AL600 was the perfect solution. Keeping my powerchair protected from the elements during transportation was very important to us. And the backlit pendant makes it just as easy to use at night!

AL600 CLASSIC PRODUCT SPECIFICATIONS

LIFTING CAPACITY	350 lb	INSTALLED WEIGHT	179 lb
COMPATIBLE VEHICLE TYPES	<ul style="list-style-type: none">- Crossovers- Minivans- SUVs- Vans	MOBILITY DEVICES	<ul style="list-style-type: none">- Manual wheelchairs- Power wheelchairs- Scooters
PLATFORM SIZE	28.5" x 38"	WHEEL CRADLE EXTENSION	28.5" x 45" (with 7" wheel cradle)
HEADROOM REQUIREMENT	32"	DEPTH REQUIREMENT	44"
SECUREMENT	4 cam buckle straps	WARRANTY	3-year
INCLUDED ACCESSORY	Wired 2-button pendant remote	AVAILABLE OPTIONS	<ul style="list-style-type: none">- 4 self-tensioning Q'Straint retractors- Battery pack- Scooter wheel cradle- Wireless key fob



Harmar.com
800.833.0478

12MAR2025 | MK-000240 A

