

UNIVERSAL INSIDE-OUT LIFT AL065

LIGHTWEIGHT LIFT FOR WORRY-FREE TRIPS



"My dad's chair looked light enough, but after lifting it in and out of his trunk so many times, I realized I was very wrong. This little lift is small and light enough for his trunk but does all the lifting we need! Couldn't be happier."

AL065 PRODUCT SPECIFICATIONS

LIFTING CAPACITY	- 200 lb	HITCH CLASS	- II or greater
INSTALLED WEIGHT	- 63 lb	ADJUSTABLE LENGTH	- 23" - 34"
ACCOMMODATES	- Manual wheelchairs - Small power wheelchairs - Scooters	POWERED	- Manual rotation - Power lifting
WARRANTY	- 3-year	SECUREMENT	- Docking device <i>(not included)</i>
AVAILABLE OPTION	- Battery pack		



Harmar.com
800.833.0478

16APRIL2024 | MKT-000213 A

