

# AL435 3-AXIS

THE ULTIMATE HOIST LIFT SOLUTION



## PREMIUM 3-AXIS INSIDE LIFT

The AL435 powered-boom maneuvers even large power wheelchairs or scooters into small vehicle openings. More versatile and flexible, for use inside a van, SUVs and pick-up trucks. The AL435 provides fully powered lifting, boom extension, and rotation for the easiest and most effortless loading. The power extend function allows for precise positioning while rotating into a tight cargo area.

- Power lifting, extension, and rotation
- Simplified base fits in a smaller footprint
- Passenger or driver side installation
- Weather-resistant hand control
- Universal base adapter included



*We wanted a lift that would fit in our minivan and be easy for us to operate when loading and unloading my husband's scooter. The AL435 was our first choice, we are able to retain use of the second row of seats, which is nice when we have our grandkids.*



**PENDANT CONTROLLER**

- Easy-to-grip base with soft, non-slip rubber coating
- Backlit buttons for night use
- Built-in LED flash light to assist with visibility while loading
- On-Off button with shut-off timer to prevent unintended activation
- Magnetic backing for convenient stowage — elevated from water in trucks
- Attractive, modern aesthetics

**AL435 3-AXIS PRODUCT SPECIFICATIONS**

<b>LIFTING CAPACITY</b>	400 lb	<b>INSTALLED WEIGHT</b>	76 lb
<b>HORIZONTAL REACH</b>	21" - 32"	<b>ADJUSTABLE HEIGHT</b>	32" - 43"
<b>WARRANTY</b>	3-year	<b>SECUREMENT</b>	Docking device <i>(not included)</i>
<b>INCLUDED</b>	<ul style="list-style-type: none"> <li>- Folding arm</li> <li>- Universal base adapter</li> </ul>	<b>ACCOMMODATES</b>	<ul style="list-style-type: none"> <li>- Manual wheelchairs</li> <li>- Power wheelchairs</li> <li>- Scooters</li> </ul>
<b>POWERED</b>	<ul style="list-style-type: none"> <li>- Power lifting</li> <li>- Power rotation</li> <li>- Power extension</li> </ul>	<b>AVAILABLE OPTION(S)</b>	<ul style="list-style-type: none"> <li>- Battery pack</li> <li>- Custom adapters</li> <li>- Custom docking devices</li> <li>- Wireless remote</li> </ul>