

# AL215 LOW PROFILE HOIST LIFT

SIMPLE, ECONOMICAL SOLUTION  
FOR LIGHTER WEIGHT MOBILITY DEVICES



Easy-to-use hand control  
with soft-touch buttons

Affordable low profile hoist lift transports and stores mobility devices inside a cargo area. The single axis lift rotates manually with powered lifting; arm height adjusts to application.

## FEATURES & BENEFITS

- Adjustable arm height
- Manual rotation
- Weather-resistant hand control
- Passenger side standard; driver side available
- Offset post allows the mobility device to be safely unloaded on the side of the vehicle

## AL215 PRODUCT SPECIFICATIONS

<b>LIFTING CAPACITY</b>	250 lb
<b>ACCOMMODATES</b>	- Manual wheelchairs - Power wheelchairs - Scooters
<b>POWERED</b>	- Manual rotation - Power lifting
<b>ARM HEIGHT</b>	26" - 34"
<b>WARRANTY</b>	3-year

<b>INSTALLED WEIGHT</b>	42 lb
<b>INCLUDED ACCESSORY</b>	Universal base adapter
<b>AVAILABLE OPTIONS</b>	- Battery pack - Custom adapters - Custom docking devices
<b>ADJUSTABLE LENGTH</b>	21" - 32"
<b>SECUREMENT</b>	Docking device <i>(not included)</i>



Harmar.com  
800.833.0478

12JULY2022 | MKT-000182 C

