

ECONOMY INSIDE LIFT AL055

SMALL BUT MIGHTY PORTABLE MICROLIFT



"When I was injured, it was a struggle getting my grandma's wheelchair in and out of her car. So, we got this little lift to help. Now that we have it, I refuse to go back to the old days of doing the lifting. I'll save that for the gym!"

AL055 PRODUCT SPECIFICATIONS

LIFTING CAPACITY	- 200 lb	INSTALLED WEIGHT	- 40 lb
ACCOMMODATES	- Manual wheelchairs - Small power wheelchairs - Scooters	ADJUSTABLE LENGTH	- 23" - 34"
AVAILABLE OPTION	- Battery pack	WARRANTY	- 3-year
POWERED	- Manual rotation - Power lifting	SECUREMENT	- Docking device <i>(not included)</i>



Harmar.com
800.833.0478

16APRIL2024 | MKT-000212 A

