

SUPERB COMFORT

EXCEPTIONAL VALUE

GENESIS
LC-150



EVOLUTION
LC-435



Sitting Pretty™ power lift recliners, excluding LC-435 XL and XXL, are FDA Class II Medical Devices designed to aid individuals with mobility impairments.



**Storage
Compartment**
Only available on
LC-435 S/M/L



**Cup
Holder**
Only available on
LC-435 S/M/L



Sitting Pretty™

New **Sitting Pretty** power lift recliners from Pride Mobility® exude style and sophistication. These three-position recliners offer superior craftsmanship and functionality. Eliminate pressure points with a chaise pad design that fully supports the legs and choose from a variety of luxurious fabrics and colors that fit your home décor. Experience at-home comfort like never before!

PRIDEMOBILITY.COM

FEATURES

- Two collections (Genesis LC-150 and Evolution LC-435), each with their own unique features and styles
- **Genesis:** Small, Medium sizes
- **Evolution:** Small, Medium, Large, XL, XXL sizes
- Easy-to-use toggle hand control
- Magnetic hand control dock on both arms (featured on all sizes except XL and XXL)
- Three-position lift and recline
- Evolution Small, Medium, and Large sizes offer an enhanced three-position experience with additional resting positions
- Chaise pad offers complete support of the legs to eliminate pressure points

OPTIONS & ACCESSORIES

- Lithium battery backup (standard on XL & XXL)
- Head and arm cover kit
- Non-skid leg levelers
- Cleaning kit

Sitting Pretty™ power lift recliners, excluding LC-435 XL and XXL, are FDA Class II Medical Devices designed to aid individuals with mobility impairments.




	GENESIS		EVOLUTION				
Model	LC-150S	LC-150M	LC-435S	LC-435M	LC-435L	LC-435XL	LC-435XXL
Size	Small	Medium	Small	Medium	Large	X-Large	XX-Large
User Height	5'3" & under	5'4" - 5'9"	5'3" & under	5'4" - 5'9"	5'9" - 6'	5'8" - 6'1"	5'8" - 6'1"
Weight Capacity	325 lbs.	325 lbs.	400 lbs.	400 lbs.	400 lbs.	500 lbs.	600 lbs.
Seat-to-Floor	18"	20"	18"	20"	22"	20"	20"
Seat Depth	18"	20"	18"	20"	22"	22"	22"
Seat Width	20"	22"	20"	22"	24"	26"	30"
Top of Back to Seat	24"	26"	24"	26"	28"	29"	29"
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Seating Type	Sinuous Spring	Sinuous Spring	Coil Pocketed Spring	Coil Pocketed Spring	Coil Pocketed Spring	Total Comfort	Total Comfort
Distance from Wall	17"	18.5"	16"	16.5"	19"	19"	19"
Overall Width	32"	33"	38.5"	40.5"	42"	39.5"	43"
Total Weight	111 lbs.	116 lbs.	136 lbs.	141 lbs.	146 lbs.	187 lbs.	196 lbs.
Warranty	2-year Limited	2-year Limited	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated
Made in America <small>of foreign and domestic parts</small>	No	No	No	No	No	Yes	Yes
Position Type	Traditional 3-position	Traditional 3-position	Enhanced 3-position	Enhanced 3-position	Enhanced 3-position	Traditional 3-position	Traditional 3-position


Fabric Options

LC-150S / LC-150M

Sophie Fabrics




Latte Sapphire



Shale Cabernet


Lexis Vinyl¹



Chestnut

LC-435S / LC-435M / LC-435L

Cosi Fabrics



Mocha Gravel Blue Diamond Syrah

Sta-Kleen¹



Chestnut


LC-435XL / LC-435XXL

Cosi Fabrics



Mocha Gravel

Sta-Kleen¹



Chestnut



pridemobility.com

800.800.8586



NOTE: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1) Additional charges will apply for this fabric.