





STIMULATE YOUR SENSES

Immerse yourself in an entertainment experience unlike any other! The Premier collection by VivaLift!® Power Recliners boasts an integrated audio system and rumble function to make you feel like you're part of your favorite movie, sporting event, or video game. Adjust the settings to your liking while enjoying premium features such as VivaLift! Tilt, heat and massage, and infinite position lounging!











100% American leather

SORRENTO COFFEE

100% American leather

SORRENTO ECLIPSE COMING SOON!







Heat and massage controls



FEATURES

- IMMERSIVE AUDIO! Speakers are built into this chair and connect to TV/Radio/Media via Bluetooth® tech.
 - Audio controls are located in the armrest.
- Heating pads in the shoulders, lumbar, and seat
- Air bladder massage system
- Hideaway wireless charger and self-leveling cup holder
- USB-A and USB-C charging ports
- Hand control features magnetic docking, memory button and tilt control
- Footrest extension
- Lithium battery backup ensures your recliner is operational when the power goes out

SPECIFICATIONS

POSITION	True-Infinite
SIZE	Medium
HEIGHT RANGE	5′5″-6′0″
WEIGHT CAPACITY	400 lbs.
SEAT-TO-FLOOR	20.5″
SEAT DEPTH	21″
SEAT WIDTH	23″
TOP OF BACK TO SEAT	26″
FULLY PADDED CHAISE	Yes
SEATING TYPE	Pocketed Coiled Springs
CHAIR HEIGHT	43″
DISTANCE FROM WALL	17"
OVERALL WIDTH	39.5″
CHAIR WEIGHT	232 lbs.
WARRANTY	7-Yr. Pro-rated

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

1. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary \pm 1". Specifications subject to change without prior notice. The information containedherein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



