



# WHILL Model F

The perfect travel chair has arrived.



# The Perfect Travel Power Chair Has Arrived.

Lightweight, foldable, maneuverable. The Model F allows you to explore your world with greater mobility and independence.

Experience the same performance and smart technology that WHILL is known for with a new foldable design.



## Foldable Design

Easily transport the Model F on buses, trains, and planes – or in the trunk of a car or taxi.\*



When not in use, the Model F can be conveniently stored in a closet or out-of-the-way corner of the room.



\*Model F is not to be used as a seat while it is being transported.

The information contained herein is correct at the time of printing, and we reserve the right to alter specifications without prior notice.

# Lightweight

Weighing only 52.9 lbs., the Model F is one of the lightest power chairs on the market today.\*\*



Top Speed



Range



Obstacle Clearance



30.7" turning radius for tight turns

# Connect Intelligently

The WHILL app allows you to check key device information and drive the Model F remotely.



Drive remotely, lock/unlock, and manage speed settings via the WHILL smartphone app.

# Drive with Peace of Mind

Stay engaged in activities and embrace your surroundings knowing the Model F will keep you safe while you ride.



Safely drives up or down 10° slopes.



\*\*Weight is 58.9 lbs. if battery is not removed. Battery weighs 6 lbs.

Note: Driving range of 12.4 miles will vary due to user weight, drive surface, terrain, ambient temperature, speed setting, and battery charge and condition.

# Customization

Choose from five colors and a great selection of accessories to match your personal style and enhance your driving experience.

## Color Options



- White   ● Black   ● Light Green
- Red   ● Light Blue



Controller can be positioned on either the left or right side.

## Accessories



Battery

Battery  
Charger

Lap Belt



Smart Key

Joystick

Basket



Cane Holder

Dust Cover

Side Bag

# Specifications

Driving Range	Up to 12.4 Miles	Maximum Speed***	3.7 MPH
Foldable	Yes	Max Incline	10°
Turning Radius	30.7"	Obstacle Clearance	1.4"
Seat Sizes Available	18"	Width Between Armrests	Adjustable: 2 positions (16", 18")
Device Width	21.9" / 23.8"	Device Length	36.8" (18.3" when folded)
Device Weight	58.9 lbs	Charging Time	~ 5 hours
App Enabled	Yes, with remote control	Air Travel Approved	Yes
Max Weight Capacity	Up to 250 lbs	Swappable Arm Covers	5 Colors
Battery	Lithium-Ion (25.3V 10Ah)	Water Resistance	IPX4

\*\*\*Acceleration and/or deceleration may change when driving on steep slopes.

## WHILL Inc.

303 Twin Dolphin Dr., 6th Floor, Redwood City, CA 94065  
 (844) 699-4455   us.info@whill.inc   www.whill.inc



Indications for Use: The intended use of the WHILL Model F powered wheelchair is to provide outdoor and indoor mobility to persons limited to a seated position that are capable of operating a powered wheelchair.

FDA Compliance: The WHILL Model F has received 510(k) clearance for marketing as a Class II medical device.

740-00226