

# MIGHTY VOYAGER



## Item# F-VGRPL

U.S. Patent  
#10,792,203

The Mighty Voyager Lift provides a whole new perspective on pool access. An ADA Compliant and UL Certified Pool Lift that is also portable allow the user to move the lift away from the side of the pool in certain cases to allow for facility events or to simply store the lift.

A robust 325 pound weight capacity is provided by separately purchasing sand that is readily available at your local pool supply or home improvement supply store. This sand ballast feature saves a considerable amount in shipping costs. The 15" opening on the ballast housing allows for easy filling and removal of the sand if necessary.

The durable stainless steel construction and powder coat finish is available in your choice of colors; standard white, gray, or tan. Seats and footrests are also available with color options in standard blue, gray, and tan.

The Mighty Voyager comes standard with a rechargeable 24v battery, wall mounted charger, and fixing anchor for ADA compliance. The Mighty Voyager is also backed by the best warranty in the industry, a 5 year structural and 5 year prorated electronics warranty.

**Shipping weight is 425 pounds and ships via truck.**

### OPTIONS/ACCESSORIES:

F-734RSA	Rescue spine board option for safe water extractions
F-MTYSLS	Sling seat option for people with greater support needs
F-MTLCB-HE	Cover, premium fade resistant, blue
F-MTLCB	Cover, standard lockable, blue
F-MTLCT	Cover, standard lockable, tan
F-MTLCG	Cover, standard lockable, gray



### ACCENT COLORS

Available in a variety of accent and powder coat colors



\*Blue seat standard (not a powder coat option)

### FEATURES:


- ADA Compliant
- UL Certified to Pool Lift Standards
- 325 lb. (147 kg) weight capacity
- 12" water draft
- ADA Fixing kit included
- Quiet & durable direct drive rotation
- Flip up armrests
- All stainless steel construction
- Durable powder coat finish in multiple colors
- Rechargeable battery & wall mount charger
- Five year pro-rated electronics warranty
- Five year structural warranty

**PROUDLY SOLD BY:**


THE MIGHTY SERIES

# MIGHTY VOYAGER

## PRODUCT SPECIFICATIONS



## APPLICABLE POOL PROFILES



\*provide measurements indicated to confirm application or use online pool lift selector

Become a member of our new Dealer and Architect page. Access to CAD blocks, CSI specs and sales materials.