

PrimeMat 2.0 Impact Reduction Fall Mats

HCPCS:




FIGURE A

- Full 72" length provides maximum coverage alongside bed
- Tough and lightweight - Mats are constructed from high quality molded Polyurethane foam inner core (210kgs/cubic meter) to provide maximum absorption and rebound to reduce impact in the event of a fall
- Passed 5 foot "Egg Drop Test" - an egg was dropped from 5 feet above the mat without damage
- Low profile with tapered edges allow equipment such as overbed tables or patient lifts to roll over easily
- Anti-slip bottom netting keeps the mat secure to help prevent accidental tripping
- Embossed non-slip PVC leather look top is not only attractive but extremely durable and puncture resistant
- Non-slip pattern helps minimize slipping and is easy to clean with mild soap and water
- Available in two colors Brown & Gray

PRODUCTS

Description

PM20BA	PrimeMat 2.0 (Brown)70x24x .75/1 EA/5 CV
--------	--