


MIGHTY VOYAGER PORTABLE LIFT

U.S. PATENT #10,792,203

*MEASUREMENT TAKEN WITH THE LIFT AT FULL LOAD



[2314mm]
91 1/8"

FOOTREST OUT

[2013mm]
79 1/4"

FOOTREST IN

[1557mm]
61 1/4"

PRODUCT WEIGHT: 315lbs [143kg]

PRODUCT WEIGHT W/ SAND: 1215lbs [551kg]


MAX WEIGHT CAPACITY: 325lbs [147kg]

MAX WATER DRAFT: *12"

SEAT HEIGHT ADJUSTMENT: 6" IN 2" INCREMENTS

VERTICAL TRAVEL DISTANCE: 48 1/2"

OUTRIGGERS
MUST BE DEPLOYED
BEFORE USE



9889 GARRYMORE LANE
MISSOULA, MT 59808
(406) 549-0769
FAX (406) 549-2602

NOTE:
*SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE



PRODUCT:		MIGHTY VOYAGER	
PART #:		F-VGRPL	
SCALE:	1:36	DATE:	11/5/2020
DRAWN BY:	KR	REVISION:	00