



GO CHAIR® MED

FDA CLASS II MEDICAL DEVICE.

(K0821) HCPCS CODE: GROUP 2: STANDARD, PORTABLE

NOW AVAILABLE IN: BLACK ONLY

- Seating: Choice of Compact, foldable 18"W x 17"D or 20"W x 18"D
 - 5 Point Swivel / 60" Lap belt / Width & Height adjustable arms
- Turning Radius: 25.5" / Overall Length: 28.5" / Overall Width: 25.5"
- Rear Wheel Drive / Limited Suspension / 1.75" Ground Clearance
- Feather touch disassembly / Heaviest Piece is 36 lbs.
- 6" Front Solid Casters / 9" Solid Drive Tires / 3" Rear Solid anti-tips
- (2) 18 Ah, 12 V Batteries / Speed Up to 3.7 mph / Range Up to 8.6 miles
 - 2 amp Off-Board Charger
- Programmable 40-amp, PG GC3 Controller

MSRP \$3,449

MAP \$2,499



FDA CLASS II MEDICAL DEVICE. PRIDE FDA CLASS II MEDICAL DEVICES ARE DESIGNED TO AID INDIVIDUALS WITH MOBILITY IMPAIRMENTS.

NOW OFFERED ON THE GO CHAIR® MED

6" FRONT CASTERS WITH EASMENT WHEELS
THAT AID IN CLIMBING OVER A 40MM
THRESHOLD



OPTIONAL SWING AWAY JOYSTICK: INLINE LEFT
OR RIGHT SIDE



THE GO CHAIR® MED IS #1 FOR PORTABILITY! WITH CONVENIENT FEATHER-TOUCH DISASSEMBLY, THE GO CHAIR MED ALLOWS YOU TO EXPERIENCE LIGHTWEIGHT TRAVEL AND INDEPENDENCE ON THE GO. COMPACT, PORTABLE AND PERFECT FOR TRAVEL.

ACCESSORIES

**WEATHER COVER
CUP HOLDER
O2 HOLDER
SAFETY FLAG
SADDLE BAG
USB CHARGER
CANE / CRUTCH HOLDER
MILITARY PATCHES**



FDA CLASS II MEDICAL DEVICE. PRIDE FDA CLASS II MEDICAL DEVICES ARE DESIGNED TO AID INDIVIDUALS WITH MOBILITY IMPAIRMENTS.

WARRANTY

- FRAME: 5-YEAR LIMITED
- ELECTRONICS: 13-MONTH LIMITED
- DRIVE MOTORS: 13-MONTH LIMITED
- BATTERY: 13-MONTH LIMITED

