

## Panther 4-Wheel Heavy Duty Scooter



HCPCS:




FIGURE A



FIGURE B



FIGURE C



FIGURE D

The Panther scooter offers modern styling and a host of excellent features resulting in a smooth and comfortable driving experience. Just some of these great features include one touch infinite tiller adjustment and delta controls for easy steering and throttle control, a plush Captains Seat that is height adjustable and features a reclining backrest with adjustable headrest, 14" alloy wheels on front and rear suspension for a smooth ride on all terrain, extra bright LED headlights and turn signals, dual rear view mirrors, and a large front basket with carry handle. For the ultimate driving experience dont settle for anything less than the Panther.

- Front and rear suspension on 14" pneumatic tires for the smoothest ride possible
- One-touch infinite tiller adjustment
- Built-in USB port for charging your mobile devices
- Extra-bright LED lighting for superb visibility
- Ergonomic delta control handles for ease of steering and throttle control
- Large front plastic basket with carry handle
- Captain's seat is height adjustable and features a fold-down, reclining backrest with adjustable headrest
- LED headlights, parking lights, and turn signals for superior visibility
- Dual rear view mirrors
- Anti-tip wheels
- Easy freewheel mode

PRODUCTS	Description
PANTHER22CS	Panther HD Scooter, 22" Capt.Seat,1/ea
SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	264 lb
Product Dimensions	56" (L) x 26" (W) x 53" (H)
Seat Dimensions	20" (D) x 22" (W)
Weight Capacity	425 lb
Top Speed	8 mph
Batteries (pair)	12V x 50AH
Battery Charger	5A Offboard
Climbing Angle	7.5°
Color	Red
Controller	Dynamic Rhino 2 120A
Ground Clearance	4"
HCPCS	N/A
Maximum Range	25 miles

**SPECIFICATIONS****Description**

Motor	24V x 800W x 4500 RPM
Seat Type	Captain - Reclining High
Turning Radius	53"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame
Wheel Size	Front: 14" x 3"   Rear: 14" x 3"