



## Phoenix HD 4 Heavy-Duty Travel Scooter

HCPCS: K0801



INSET A




INSET B



INSET C



INSET D



INSET E



INSET F



INSET G

- Quick and easy disassembly
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- Quick connect batterie
- Ergonomic throttle control
- Height-adjustable swivel seat with two-tone upholstery
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Headlight
- Swivel seat with fold-down backrest
- Easy free wheel operation

PRODUCTS	Description
PHOENIXHD4	Phoenix HD 4Whl Travel POV 17.5 Seat, 1/cs
PHOENIXHD4-20	Phoenix HD 4Whl Travel POV 20 Seat, 1/cs

SPECIFICATIONS	Description
HCPCS	K0801
Weight Capacity	350 lbs.
Top Speed	4 mph
Maximum Range	15 miles
Turning Radius	54"
Climbing Angle	6°
Ground Clearance	3.75"
Dimensions	41.5" (L) x 22.5" (W)
Floor-to-Seat Height	20" - 22"
Seat Type	Folding
Seat Dimensions	17.5" (W) x 17" (D)
Motor	24V x 350W x 3600 rpm
Batteries (pair)	12V x 20AH
Charger	1.5A Offboard
Brakes	Electromagnetic
Front Wheels	9" x 3" Flat-free
Rear Wheels	9" x 3" Flat-free
Base Weight	77 lbs.
Battery Weight	31 lbs.

SPECIFICATIONS	Description
Seat Weight	21 lbs.
Total Weight	129 lbs.
Heaviest Piece	46 lbs. (Front Section)
Warranty on Frame	Lifetime
Warranty on Electronic Controller	24 Months
Warranty on Drive Train Components	24 Months
Wheel Base Dimensions (L x W)	30.7" x 19.1"
Warranty on Batteries	12 Months
Wheel Base Dimensions (L x W)	30.7" x 19.1"

ACCESSORIES	Description
SF8010	Oxygen Tank Holder, 1/ea
SF8020	Rear Basket, 1/ea
SF8030	Cane/Crutch Holder, 1/ea
SF8050	Sun Shade, 1/ea