

## Enhance Mobility



# TRIAXE

## SPORT



The peppy TRIAXE Sport Scooter can travel at speeds up to 12mph and distance to 35 miles on a single charge. This light weight travel scooter folds flat for storage and transport.

## TOUR



This lightweight, folding travel scooter features rear wheel drive, paddle acceleration and electromagnetic brakes. The aluminum frame adds strength, reduces weight and eliminates corrosion to provide years of reliability.

<b>Construction</b>	Aluminum Frame	Aluminum Frame
<b>Folds Flat for Storage</b>	Yes	Yes
<b>Weight (without Battery)</b>	44lbs	58 lbs
<b>Speed</b>	Up to 12mph	Up to 6mph
<b>Weight Capacity</b>	350lbs	375lbs
<b>Battery</b>	48V 12AH (travel battery available - 48V 6AH)	24V 12amp (12AH)
<b>Brake</b>	Dual Rear Disc Brakes	Electric Brake System
<b>Length</b>	34"	38"
<b>Width</b>	20 ½"	19"
<b>Turning Radius</b>	33"	27"
<b>Range</b>	Up to 35 miles	Up to 18 miles
<b>Drive System</b>	Front Wheel Drive Brushless Motor	Rear Wheel Drive Transaxle with Free Wheel Mode
<b>Colors Available</b>	Metallic Red, Light Metallic Blue, Green, Black, Silver, Purple (available soon)	Fiery Red, Bright Blue
<b>Seat</b>	13" standard, 18" option	13" standard, 18" option
<b>Tires</b>	8" Solid Front 9.5 Pneumatic Rear Tires with Anti-Flat	8" Solid Front 9.5 Pneumatic Rear Tires with Anti-Flat

Enhance Mobility

[sales@enhancemobility.com](mailto:sales@enhancemobility.com)

800-296-4151