

ZooMe Auto-Flex Folding Travel Scooter



INSET A



INSET B



INSET C



INSET D



INSET E



INSET G



INSET H



INSET I

- Fully electronic folding and unfolding in 15 seconds or less with a wireless key fob
- Aluminum alloy frame is lightweight for superior portability
- Airline safe lithium battery
- High-gloss shrouds for a protective finish
- Easy throttle control
- Delta tiller for users with limited dexterity
- Auto-folding seat with padded, adjustable armrests
- Flat-free, non-marking tires for worry-free travel
- Easy-to-adjust tiller places the scooter's controls in the ideal position
- Secondary manual-release levers and toggle switch for folding and unfolding
- Front and rear reflectors
- 2 ignition keys and 1 key fob included
- Anti-tip wheels for additional safety

PRODUCTS

PRODUCTS	Description
FLEX-AUTO	Electric Blue, 1/ea

SPECIFICATIONS

SPECIFICATIONS	Description
Weight Capacity	300 lbs.
Top Speed	4 mph
Estimated Range	13 miles
Climbing Angle	6°
Turning Radius	47.2"
Ground Clearance	1.6"
Overall Dimensions	37" (L) x 17" (W) x 32.5" (H)
Overall Dimensions Folded	17" (L) x 17" (W) x 23.5" (H)
Seat Dimensions	16" (W) x 14" (D)
Motor	24V x 120W
Controller	PG S Drive 45A
Batteries	25.9V x 10AH Lithium
Battery Charger	2A Offboard
Brakes	Electromagnetic
Freewheel Mode	Yes
Front Wheels	6" x 1.5"
Rear Wheels	7" x 2.35"
Anti-Tip Wheels	2" x .5"
Base Weight	55 lbs.
Battery Weight (Pair)	5 lbs.

SPECIFICATIONS

Description

Heaviest Piece Weight	55 lbs.
Total Weight	60 lbs.
Warranty on Frame	Limited Lifetime
Warranty Motor Electronics Battery	24 24 12 Months