

## ZooMe Flex Folding Travel Scooter Travel Scooter



INSET A



INSET B



INSET C



INSET D



INSET E



INSET F

- Ultra-lightweight for superior portability
- Quick and Easy "Flex" folding design for transport and storage in one easy motion
- High-gloss shrouds for a protective finish
- Easy throttle control
- Delta tiller for users with limited dexterity
- Auto-folding seat with padded, adjustable armrests and fold down back
- Flat-free, non-marking tires for worry-free travel
- Easy-to-adjust tiller angle places the scooter's controls in the ideal position
- Anti-tip wheels for additional safety

### PRODUCTS

PRODUCTS	Description
FLEX	Electric Blue, 1/ea

### SPECIFICATIONS

SPECIFICATIONS	Description
Weight Capacity	275 lbs.
Top Speed	4 mph
Estimated Range	13 miles
Turning Radius	53"
Climbing Angle	6°
Ground Clearance	2"
Overall Dimensions	37" (L) x 18" (W) x 36" (H)
Seat Dimensions	16.75" (W) x 15" (D)
Motor	24V x 120W
Controller	PG S Drive 45A
Batteries	25.9V x 10AH Lithium
Battery Charger	2A Offboard
Brakes	Electromagnetic
Freewheel Mode	No
Front Wheels	6" x 1.5"
Rear Wheels	7" x 2.35"
Anti-Tip Wheels	2" x .5"
Base Weight	53.6 lbs.
Battery Weight (Pair)	5 lbs.
Heaviest Peice Weight	53.6 lbs
Total Weight	58.6 lbs.
Warranty	24 months on Electronics, 12 months on Batteries, Limited Lifetime on Frame

### PARTS

PARTS	Description
FLEX-01	UPPER PANEL ASSY, 1/EA

PARTS	Description
FLEX-02	THROTTLE ASSY, 1/EA
FLEX-03	GAS SPRING ASSY, 1/EA
FLEX-04	HEADLIGHT ASSY, 1/EA
FLEX-05	THROTTLE POTENTIOMETER, 1/EA
FLEX-06	QUICK RELEASE CLAMP, 1/EA
FLEX-07	FRONT WHEEL, 1/EA
FLEX-08	STEERING KNUCKLE, LEFT, 1/EA
FLEX-09	STEERING KNUCKLE, RIGHT, 1/EA
FLEX-10	HUB CAP, 1/EA
FLEX-11	WHEEL/STEM BEARING, 1/EA
FLEX-12	FRONT COVER, 1/EA
FLEX-13	TIE ROD WITH ENDS, 1/EA
FLEX-14	HANDLE CONNECTOR, 1/EA
FLEX-15	STEERING CONNECTOR, 1/EA
FLEX-16	BATTERY CASE EMPTY, 1/EA
FLEX-17	FRONT PEDAL COVER, 1/EA
FLEX-18	ROTATING GEAR ASST, FRONT, 1/EA
FLEX-20	RELEASE CABLE, 1/EA
FLEX-19	SUPPORT ROD, FRONT, 1/EA
FLEX-21	RELEASE LEVER, 1/EA
FLEX-22	BATTERY, LITHIUM, 1/EA
FLEX-23	SEAT SUPPORT ASSY, LEFT, 1/EA
FLEX-24	SUPPORT ROD ASSY, REAR, 1/EA
FLEX-25	REAR WHEEL, 1/EA
FLEX-26	REAR PEDAL COVER, 1/EA
FLEX-27	BATTERY CHARGER, 1/EA
FLEX-28	POWER HARNESS, 1/EA
FLEX-29	CONTROLLER, 1/EA
FLEX-30	ROTATING GEAR ASSY, REAR, 1/EA
FLEX-31	BATTERY TERMINAL (4 EA), 4/EA
FLEX-32	ANTI TIP WHEEL, 1/EA
FLEX-33	CHARGER HARNESS, 1/EA
FLEX-34	REAR COVER, 1/EA
FLEX-35	MOTOR, TRANS, BRAKE, 1/EA
FLEX-36	SEAT SUPPORT ASSY, RIGHT, 1/EA
FLEX-37	ARMPAD W/SCREWS, 1/EA
FLEX-38	BACK UPHOLSTERY, 1/EA
FLEX-39	BACK CUSHION, 1/EA
FLEX-40	SEAT UPHOLSTERY, 1/EA
FLEX-41	SEAT CUSHION, 1/EA
FLEX-42	MAIN WIRING HARNESS, 1/EA
FLEX-43	Left Hinge Cover, 1/EA
FLEX-44	Right Hinge Cover, 1/EA
FLEX-45	Armrest Locking Knob, 1/EA
FLEX-46	Speed Controller, 1/EA
FLEX-47	On/Off Switch, 1/EA
LDR410115-03	Folding Locking Knob, 1/EA
FLEX-48	Armrest Assembly, 1/ea
FLEX-49	Armrest Assembly, 1/ea
FLEX-50	Docking Station, 1/ea